

RV-SH-320



REVOLVE

REVOLVE DIGITAL HOTPLATE STIRRER



PRODUCT SPECIFICATION

Motor	DC
No of stirring position	1
Max stirring quantity	10L
Speed Range (RPM)	200-2200RPM
Speed control	10RPM/Step
Heating Temp range	Room temp to 320c
Setting temp	5-320c
Safe temp	345
Temperature control	1c/step
Time range	1-999 min & infinite
Recommended stirrer bar length	25mm
weight	2kg

GET IN TOUCH



Revolve Scientific Equipment Trading LLC
Abu Dhabi, UAE, info@revolvescientific.com
sales@revolvescientific.com, service@revolvescientific.com
+971 54 710 8986, +971 54 710 8987, +971 54 710 8988

WWW.REVOLVESCIENTIFIC.COM