


## ACCESSORIES

revolvescientific.com

## MICRO CENTRIFUGE WITH ADJUSTABLE SPEED AND TIMER


REVOLVE SCIENTIFIC

# RV- MICRO-06




### TUBE REDUCTION ADAPTERS

For 0.2/0.5ml tubes



### PCR ADAPTER

For 02 x 8 x 0.2 ml



Revolve microfuge has reliable performance, intelligently designed chassis with user friendly interface

revolvescientific.com

## ASK FOR DEMO

Revolve Microfuge centrifuge comes equipped with digital display with adjustable speed and timer. Equipped with high quality suction pads, Revolve Microfuge has quite operation, can be used for extended runs at full speed with stable foot print. The fresh, red color of Revolve Microfuge micro centrifuge as well as the elegant Danish design will bring life and contrast to your lab environment. Revolve Microfuge has the centrifugation power of 6000 rpm / 2000 g. Revolve Microfuge are suitable for 1.5 and 2.0 mL tubes. Smaller micro centrifuge tubes can be inserted using the accessory adaptors.

REVOLVE



# Specification

**Rotor Capacity** 8 Place Closed Rotor for 1.5/2.0 ml microtubes

**Maximum RPM / RCF (g)** 6000 / 2000 g

**Display** Digital

**Accuracy**  $\pm 50$  RPM

**Noise Level** Yes

**Variable Speed** 2000 – 6000 RPM in steps of 100 RPM

**Timer** 1 to 25 mins & infinite

**Quick Spin** Yes, Lid activated On/Off

**Dimensions (WxDxH)** 162 x 157 x 115 mm

**Rotors and Adaptors (included in standard pack)**

8 Slot x 1.5/2.0 ml Closed Rotor

PCR Strip Rotor (2 x 8 x 0.2 ml)

Reduction adaptors for 0.2/0.5 ml Microtubes

•

## ORDERING INFORMATION

- RV-MICRO-06 Revolve Microfuge with adjustable speed and timer
- RV-MICRO-A8M 8 x 1.5 / 2 ml tubes slot magnetic rotor
- RV-MICRO-B2M PCR strip magnetic rotor
- RV-MICRO-D1 Reduction Adaptors of 0.2 ml tubes
- RV-MICRO-D2 Reduction Adaptors of 0.5 ml tubes

DO  
YOU  
KNOW

REVOLVE MICROFUGE HAS  
VARIABLE TIMER SETTING  
FROM 1 MIN TO 25 MIN AND  
INFINITE MODE

**REVOLVE**

## WEBSITE

[WWW.REVOLVESCIENTIFIC.COM](http://WWW.REVOLVESCIENTIFIC.COM)



**REVOLVE**

[www.revolvescientific.com](http://www.revolvescientific.com)