

## FEATURES

- Motorless magnetic stirrer with very slim profile design
- No moving parts – Truly maintenance free
- Adjustable speed from 15 to 1500 RPM in steps of 50 RPM
- Timer range 1- 99 mins & infinite mode
- Stirring Capacity of upto 800 ml
- IP 65 compliant design (Water spill resistant)
- Flat surface – Easy to wipe and clean
- High strength PTFE Magnetic bar (25 x 8 mm)

## ASK FOR DEMO

Revolve Magnetic Stirrer is a compact magnetic stir plate with adjustable speed up to 1500 rpm. The top surface of Revolve Magnetic Stirrer is made of robust stainless steel and is resistant towards moderate contamination of harsh chemicals. Revolve Magnetic Stirrer is also water spill resistant (IP 65 compliant). Equipped with a modern real-time digital display, magnetic stir plate shows the current speed and time settings. The touch control pad of magnetic stirrer enables an easy and quick operation and the flat surface makes cleaning extremely easy. Thanks to the last run memory feature, you can effortlessly continue your work with Revolve Magnetic Stirrer.

Revolve Magnetic Stirrer plate has the additional function of pulse mode, improving the mixing process, by stirring for 30 seconds in each direction. Only 1.65cm thick, Revolve Magnetic Stirrer is one of the slimmest profile stirrers on the market. With its 600 grams, Revolve Magnetic Stirrer is also one of the lightest. Revolve Magnetic Stirrer plate is maintenance free. The total load of Revolve Magnetic Stirrer is 800mL. The best stirring effect can be obtained with a bread base 250ml beaker. Revolve Magnetic Stirrer starts mixing at the speed of 15 rpm and can go up to 1500 rpm with speed increments of 50 rpm.

revolvescientific.com

## DIGITAL MAGNETIC STIRRER



REVOLVE SCIENTIFIC

# RV-MS-15

Chemical Resistant SS top  
Adjustable Speed and Time

revolvescientific.com



# REVOLVE

DO YOU KNOW



## Specification

- Maximum speed 1500 rpm
- Speed setting 15 to 1500 rpm (increments of 50 rpm)
- Load capacity 800 mL
- Recommended stir bar length 25 mm
- Dimensions 124 x 191 x 16.5 mm
- Weight 600 grams
- Pulse function YES
- Timer function 1 to 99 minutes and infinite mode



## ORDERING INFORMATION

- RV-MS-15 Revolve Magnetic Stirrer 1500 RPM Adjustable Speed and Time

Revolve Magnetic Stirrer has Pulse Mode



## WEBSITE

[WWW.REVOLVESCIENTIFIC.COM](http://WWW.REVOLVESCIENTIFIC.COM)

**REVOLVE**

**REVOLVE**  
[www.revolvescientific.com](http://www.revolvescientific.com)

- PROGRAMMABLE PULSE MODE WITH RUN PAUSE RUN FEATURE
- HIGHLY TEMPERATURE STABLE & CHEMICALLY INERT COMPOSITE FIBRE BODY
- AUTO START POWER INTERMITTENT MODE
- HOT SURFACE WARNING & ERROR INDICATION FEATURES
- 6 INTELLIGENT HEATING MODE
- KEYPAD LOCK TO PREVENT ACCIDENTAL CHANGES OF PARAMETERS