

FEATURES

revolvescientific.com

REVOLVE PLATFORM SHAKER

REVOLVE SCIENTIFIC
RV-PS-20

- *Microprocessor Controlled*
- *50 to 300 RPM*
- *Timer Setting*
- *Load Capacity 7.5Kg*
- *Orbital & Linear Reciprocating Motion*
- *2 Years of Warranty*



1 to 99 hours and 59 minutes, and the speed can be programmed from 50 to 300 RPM.



ASK FOR DEMO

Revolve Platform Shaker is universally compatible with almost all widely used vessels and comes with several attachment platforms.

Power failure prevention is built into the Revolve Platform Shaker, and the choice of Auto Power On will guarantee continued functioning for the remaining period following restart.

REVOLVE




**DO
YOU
KNOW**

Set motion in either clockwise & counter clockwise direction

Specification

- Motor Type
Brushless DC Motor
- Time Setting Range
Upto 99:59 hrs
- Platform Size
330 x 330 mm
- Programmable Feature
6 programs with each program having 9 sub-programs
- Program Looping
Yes
- Dimensions of Equipment
408 x 350 x 99 mm
- Maximum weight capacity
7.5 Kg



ORDERING INFORMATION

- RV-PS-20 Revolve Platform Shaker
- RV-PS-U Revolve Universal Platform with adjustable bars (Flask capacity 9 x 250mL or 6 x 500mL)
- RV-PS-C9-250 Revolve Base plate with 9 Stainless steel clamps : 9 x 250 mL
- RV-PS-C9-500 Revolve Base Plate with 9 stainless steel clamps for 9 x 500mL flasks



- MICROPROCESSOR CONTROLLED FOR EFFICIENT PERFORMANCE EVEN UNDER HIGH LOAD
- VARIABLE SPEED SELECTION UPTO 300 RPM (PS 20)
- IP21 COMPLIANT FOR MAXIMUM SAFETY AND DURABILITY
- UNIQUE LOOP FUNCTION
- SET MOTION IN EITHER CLOCKWISE & COUNTER CLOCKWISE DIRECTION

WEBSITE

WWW.REVOLVESCIENTIFIC.COM

REVOLVE

REVOLVE

www.revolvescientific.com