

Revolve Vortex Mixer 4500rpm.
Microprocessor Controlled.
Variable Speed up to 4500 rpm.
Chemical Resistant.
2 Years Warranty.
Can be operated from 4-40°C ideal
for cold rooms and incubators..



REVOLVE



**MINI VORTEX MIXER
4500 RPM**

REVOLVE

REVOLVE Mini Vortex Mixer comes with a powerful brushless DC motor and delivers speed up to 4500 RPM. With an orbital motion of 4.5 mm it provides thorough mixing and it can accommodate plastic tubes & vials up to a 30 mm diameter.

RV-VM-45F



Motion	Orbital
Speed Range (RPM)	0 – 4500
Orbital Diameter	4 .5 mm
Operating Mode	Touch
Max tube diameter	30 mm
Motor	Brush-less DC motor
Max. Load Capacity	500 gm
Permissible ambient temperature	5 - 40 °C
Permissible relative moisture	80%

GET IN TOUCH

Revolve Scientific Equipment Trading LLC
Abu Dhabi, UAE, info@revolvescientific.com
sales@revolvescientific.com, service@revolvescientific.com
+971 54 710 8986, +971 54 710 8987, +971 54 710 8988
WWW.REVOLVESCIENTIFIC.COM