

## FEATURES

- Max capacity 4 x 100 ml (swing out & fixed rotor)
- Maximum tube capacity 32 x 10 ml (swing out)
- Microprocessor controlled centrifuge with variable speed and Time setting and also with last run memory function
- Adjustable ramp & brake speeds from (up to 9 speed selection)
- **Automatic Imbalance Detection**
- **2 Years Warranty**

revolvescientific.com

REVOLVE UNIVERSAL CENTRIFUGE WITH 15000 RPM



REVOLVE SCIENTIFIC  
**RV-UC-15**

15000 RPM / 22388 g (RPM and RCF varies as per rotor)

## ASK FOR DEMO

Revolve Universal Centrifuge RV-UC-15 is equipped with industry standard features and more. Revolve Universal Centrifuge has max speed of 15000 RPM/ 22388RCF and max load capacity of 4 x 100 ml. Revolve Universal Centrifuge can be used with variety of fixed angle and swingout rotors. Revolve Universal Centrifuge has modern auto rotor detection, imbalance detection and Lid/Lock safety Features.

revolvescientific.com



**REVOLVE**

**DO YOU KNOW**

Revolve Universal Centrifuge has short spin function

- ADJUSTABLE RAMP & BRAKE SPEEDS FROM (UP TO 9 SPEED SELECTION)
- IMBALANCE DETECTION AND LID LOCK SAFETY FEATURE
- USER-PROGRAMMABLE (USER CAN SET UP TO 99 PROGRAMS WITH 4 LINE PROGRAMS)
- BRUSHLESS DC MOTOR FOR MAINTENANCE-FREE LONG LIFE

**REVOLVE**

## Specification

Rotor Capacity (Max) 4 x 100 ml  
32 x 10 ml (Swing out) & 30 x 15 ml (fixed rotor)  
Max RPM / RCF (xg) 15000 RPM / 22388 g (RPM and RCF varies as per rotor)  
Motor type Brushless DC  
Display Digital display  
Timer Setting 30 sec to 999 mins & infinite mode  
Lid Lock Safety Yes  
Imbalance Detection Yes  
Adjustable Ramp & Brake Yes (selection from 1 to 9)  
Weight 23 kg (without rotor)  
Rotors (not included) Multiple rotor options (fixed angle and swing out)

For additional Rotor options c refer to our website

**WEBSITE**

[WWW.REVOLVESCIENTIFIC.COM](http://WWW.REVOLVESCIENTIFIC.COM)

## ORDERING INFORMATION

- RV-UC-15 Revolve Universal Centrifuge with 15000 RPM (without rotor)
- RV-PCR-F-4X8 Revolve 4X8X0.2 PCR rotor
- RV-2-F-48 Revolve Fixed angle rotor 44x2ML
- RV-100-F-4 Revolve 4 x 100 mL fixed angle rotor, RCF 2559 x g
- RV-50-F-4 Revolve 4 x 50 mL fixed angle rotor, RCF 2644 x g
- RV-15-F-16 Revolve 16 x 15 mL fixed angle rotor, RCF 3144 x g
- RV-100-S-4 Revolve 4 x 100 mL swing rotor with Safety Cap, RCF 3164 x g



**REVOLVE**

[www.revolvescientific.com](http://www.revolvescientific.com)