



SAFE



INTELLIGENT



PRECISION




**REVOLVE**

**HEATING MANTLE WITH STIRRER**



The sturdy elegant HMC series of precision heating mantles by Revolve are temperature controlled that precisely control fluid temperature upto an accuracy of  $\pm 1$  c. They have automatic power cutoff option after exceeding set temperature.



Durable and elegant HMS series of precision heating mantles are:

- Temperature controlled that precisely control fluid temperature up to an accuracy of  $\pm 1$  c.
- Safe temperature settings to prevent accidents.
- It is designed to accommodate all types of flasks.
- Flexible yarn heating element.
- LED display indicating the process status.
- Insulated, grounded metal body with a protective coating
- Corrosion-resistant stainless steel top
- Pre-threaded hole at the back for attaching accessories

HEATING MANTLE WITH STIRRER



Model	RV-HMS-500	RV-HMS-1000
Flask Capacity	500 ml	1000 ml
Temperature control	PID Control and digital energy regulator	PID Control and digital energy regulator
Unit External Dimensions (WxDxH)	200mm x 200mm x250mm	265mm x 235mm x 300mm
Working Temperature range	Ambient to +5°C to 400°C	Ambient to +5°C to 400°C
Display	7 Segment LED	7 Segment LED
Temperature Accuracy	±1.0°C	±1.0°C
Timer	1-999 Minutes	1-999 Minutes
External probe	PT100	PT100
Speed range	300-1500 Rpm	300-1500 Rpm

GET IN TOUCH



Revolve Scientific Equipment Trading LLC  
Abu Dhabi, UAE, info@revolvescientific.com  
sales@revolvescientific.com, service@revolvescientific.com  
+971 54 710 8986, +971 54 710 8987, +971 54 710 8988

WWW.REVOLVESCIENTIFIC.COM