

CHILLER TANK CAPACITY 10L



REVOLVE

REVOLVE

The Revolve chiller features a temperature range from -10°C to ambient -5°C and boasts an ergonomic, compact design that allows it to be positioned at eye-level on the desk. It includes a large, bright LED display with reliable microprocessor PID temperature control, ensures economic and environmentally friendly operation, and offers a simple operation with an easy-to-navigate menu.





PRODUCT SPECIFICATION

Tank capacity (L)	10
Tank size (mm)	325mm x 265mm x 200mm
Cooling capacity (Watt)	1000
Weight	52
Pump flow rate (LPM)	12
Display	LED
Temperature stability	±1°C
Temperature range	-10 to 25°C
Flow rate	12 LPM
Overall size (LxWxH) mm	594mm x 396mm x 723mm
Ordering information	RV-CHL-10

COOLING CAPACITY CHART

TEMP (°C)	20	10	0	-10
COOLING CAPACITY (W)	1000	800	600	350

GET IN TOUCH



Revolve Scientific Equipment Trading LLC
Abu Dhabi, UAE, info@revolvescientific.com
sales@revolvescientific.com, service@revolvescientific.com
+971 54 710 8986, +971 54 710 8987, +971 54 710 8988

WWW.REVOLVESCIENTIFIC.COM