



REVOLVE

RV-AB-224
REVOLVE ANALYTICAL BALANCE



REVOLVE ANALYTICAL BALANCE

High-precision weighing for laboratory, pharmaceutical, or industrial applications

Features include:

- Capacity 0-220g/0.0001g
- Platform Size 90 mm
- Stable Time 1.5s
- Preheating 30-60 mins
- Weighing Units: g/oz/ct/mg
- Min Weight 10 mg
- Repeatability $\pm 0.1\text{mg}$
- Linearity $\pm 0.2\text{mg}$
- RS232 and USB bidirectional interface
- Electromagnetic Force Sensor
- Super Bright Detachable Displayer
- Multi-Unit Weighing Function, Customizable Units
- Glass Windshield-Antistatic Material
- Fast Sensitivity Setting
- Below Weighing System
- Time Display
- The Weighing Results can be Exported to Excel
- Built in Stand Alone Dual Program (0.1mg)
- Power AC 110V-240V
- Weighing Function: Basic weighing, Dynamic weighing, Piece Counting, Percent Weighing, Density Determination, Pipette Calibration Function
- GMP/GLP (Touch Screen)

Improved user experience with upgraded premium features, including:

- Quick access to menu options via a large 7-inch high-resolution full-touch display with a sleek, modern interface.
- Multiple convenient applications to enhance weighing efficiency.
- Smooth connectivity options for easy data transfer, interface compatibility, and GLP/GMP-compliant printouts.





GET IN TOUCH

REVOLVE

Revolve Scientific Equipment Trading LLC
Abu Dhabi, UAE, info@revolvescientific.com
sales@revolvescientific.com, service@revolvescientific.com
+971 54 710 8986, +971 54 710 8987, +971 54 710 8988

WWW.REVOLVESCIENTIFIC.COM