

REVOLVE



RV-AB-224I  
REVOLVE ANALYTICAL BALANCE  
WITH INTERNAL CALIBRATION



## REVOLVE ANALYTICAL BALANCE WITH INTERNAL CALIBRATION

Revolve analytical balances with internal calibration offer high precision, automatic calibration, and user-friendly operation, making them ideal for laboratory applications.

### Features include:

- Capacity 0-220g/0.0001g
- Platform Size 90 mm
- Stable Time 1.5s
- Preheating 30-60 mins
- Weighing Units: g/oz/ct/mg
- Min Weight 10 mg
- Repeatability  $\pm 0.1\text{mg}$
- Linearity  $\pm 0.2\text{mg}$
- RS232 and USB bidirectional interface
- Electromagnetic Force Sensor
- Super Bright Detachable Displayer
- Multi-Unit Weighing Function, Customizable Units
- Glass Windshield-Antistatic Material
- Fast Sensitivity Setting
- Below Weighing System
- Time Display
- The Weighing Results can be Exported to Excel
- Built in Auto Calibration System
- Built in Stand Alone Dual Program (0.1mg)
- Power AC 110V-240V
- Weighing Function: Basic weighing, Dynamic weighing, Piece Counting, Percent Weighing, Density Determination, Pipette Calibration Function
- GMP/GLP (Touch Screen)

### Improved user experience with upgraded premium features, including:

- Quick access to menu options via a large 7-inch high-resolution full-touch display with a sleek, modern interface.
- Multiple convenient applications to enhance weighing efficiency.
- Smooth connectivity options for easy data transfer, interface compatibility, and GLP/GMP-compliant printouts.





## GET IN TOUCH

# REVOLVE

Revolve Scientific Equipment Trading LLC  
Abu Dhabi, UAE, [info@revolvescientific.com](mailto:info@revolvescientific.com)  
[sales@revolvescientific.com](mailto:sales@revolvescientific.com), [service@revolvescientific.com](mailto:service@revolvescientific.com)  
+971 54 710 8986, +971 54 710 8987, +971 54 710 8988

[WWW.REVOLVESCIENTIFIC.COM](http://WWW.REVOLVESCIENTIFIC.COM)