

PRECISION ELECTRONIC BALANCE
High Precision Strain Gauge Load cell
Capacity 220g/0.01g
High Contrast LCD Display
RS232C interface
Stable Time 1.5s
Units g/oz/ct/mz
Automatic Fault Detection
Overload Protection
TARE function
Unit Conversion, Counting, Percentage etc.
Platform Size 170 x 170 mm



REVOLVE



RV-PB-202E
REVOLVE PRECISION
ELECTRONIC BALANCE

REVOLVE



GET IN TOUCH

Revolve Scientific Equipment Trading LLC
Abu Dhabi, UAE, info@revolvescientific.com
sales@revolvescientific.com, service@revolvescientific.com
+971 54 710 8986, +971 54 710 8987, +971 54 710 8988
WWW.REVOLVESCIENTIFIC.COM