

REVOLVE



RV-PB-203E
REVOLVE PRECISION BALANCE,
READABILITY: 0.001-200G



REVOLVE PRECISION BALANCE, READABILITY: 0.001-200G

Revolve Precision Electronic Balance has high-precision strain gauge loadcell which delivering excellent cost vs performance ratio, Revolve Precision Electronic Balance is equipped with High contrast LCD display for easy reading of the parameters. Revolve Precision Electronic Balance has inbuilt automatic fault detection, overload protection giving user with peace of mind operation. With Revolve Precision Balance you can do Unit conversion, counting, percentage, Tare function etc.

Main Features:

- High contrast LCD display
- Big size stainless steel scale pan
- Glass windshield
- Automatic fault detection, overload protection
- Unit conversion, counting, percentage, etc.
- Tare function
- RS232 interface (Option USB interface)
- Below weighing facility (Optional)
- Dual display (Optional)
- 4 x AA battery (Optional)
- Capacity: 0-200g
- Readability: .001g
- Platform Size 90 mm
- Repeatability ($\leq$): $\pm 2\text{mg}$
- Linearity ($\leq$): $\pm 2\text{mg}$
- RS232 interface
- Power: 110V/220-240V AC





GET IN TOUCH

REVOLVE

Revolve Scientific Equipment Trading LLC
Abu Dhabi, UAE, info@revolvescientific.com
sales@revolvescientific.com, service@revolvescientific.com
+971 54 710 8986, +971 54 710 8987, +971 54 710 8988

WWW.REVOLVESCIENTIFIC.COM